

**Applied Physics****Lab Report#3****SUBMITTED BY:**

NAME	REG ID
Muhammad Maaz Abbasi	242308
Ahmed Owais	243394
Muhammad Haseeb	242292
Ahmad Saleem	242348
Syed Wasif Ali Shah	242334
Syed Ahtsham Hussain Shah	242364

**SUBMITTED TO:****Mr. Nazir Ahmed Malik**

## **EXP. 1**

### **Color Addition**

---

#### **Objective**

To understand and demonstrate the concept of additive color mixing using primary colors of light (Red, Green, and Blue).

#### **Apparatus**

- Three primary color LED light sources (Red, Green, Blue)
- Projection screen or white surface
- Power supply and controllers for LEDs
- Ruler or measuring scale

#### **Procedure**

1. Set up the three LED sources at equal angles directed towards a white screen.
2. Turn on one LED at a time to observe the individual color.
3. Switch on two LEDs at once to observe the resultant color due to mixing.
4. Finally, turn on all three LEDs simultaneously and observe the white light formed.

#### **Observations and Calculations**

LEDs On	Resultant Color
---------	-----------------

Red + Green	Yellow
-------------	--------

Red + Blue	Magenta
------------	---------

Green + Blue	Cyan
--------------	------

Red + Green + Blue	White
--------------------	-------

**Conclusion**

The experiment verifies the principle of color addition. Mixing two primary colors gives a secondary color, and combining all three results in white light.

---

## EXP. 2

### Prism

**Objective**

To study the refraction and dispersion of light through a glass prism.

**Apparatus**

- Glass prism
- White light source
- Screen
- Protractor
- Ruler

**Procedure**

1. Place the prism on a flat surface and direct the white light towards one face.
2. Mark the incident and refracted rays.
3. Measure the angle of deviation.
4. Observe the dispersion spectrum on the screen.

**Observations**

- White light splits into a spectrum of seven colors.
- The angle of deviation increases with wavelength.

**Conclusion**

The experiment demonstrates that prisms disperse light due to different refractive indices for each wavelength, validating the phenomenon of dispersion.

---

## **EXP. 3 (Reflection)**

**Objective**

To verify the laws of reflection.

**Apparatus**

- Plane mirror
- Ray box or laser pointer
- Protractor
- White paper
- Pencil

**Procedure**

1. Fix the mirror vertically on the paper.
2. Direct an incident ray at various angles.
3. Trace the incident and reflected rays.
4. Measure the angles using a protractor.

**Observations**

**Angle of Incidence ( $\theta_i$ )    Angle of Reflection ( $\theta_r$ )**

20°

20°

**Angle of Incidence ( $\theta_i$ )   Angle of Reflection ( $\theta_r$ )**

40°

40°

60°

60°

**Conclusion**

The angle of incidence is equal to the angle of reflection, validating the law of reflection.

---

## EXP. 4

### Snell's Law

**Objective**

To verify Snell's Law of Refraction.

**Apparatus**

- Glass slab or semicircular glass block
- Ray box
- Protractor
- White paper

**Procedure**

1. Shine a light beam at various angles into the glass.
2. Trace incident and refracted rays.
3. Measure the angles and compute sine values.
4. Plot  $\sin(\theta_i)$  vs  $\sin(\theta_r)$ .

**Observations and Calculations**

$\theta_i$	$\theta_r$	$\sin(\theta_i)$	$\sin(\theta_r)$
------------	------------	------------------	------------------

20°	13°	0.342	0.225
-----	-----	-------	-------

$\theta_i$   $\theta_r$   $\sin(\theta_i)$   $\sin(\theta_r)$

40° 25° 0.643 0.423

**Conclusion**

The plot is linear, confirming Snell's Law:  $\sin(\theta_i)/\sin(\theta_r) = \text{constant}$  (refractive index).

---

## EXP. 5

### Total Internal Reflection

**Objective**

To observe the conditions for Total Internal Reflection (TIR).

**Apparatus**

- Semicircular glass block
- Ray box
- Protractor

**Procedure**

1. Direct a light ray from the flat face into the curved surface.
2. Increase the angle of incidence gradually.
3. Identify the critical angle where the refracted ray disappears.

**Observations**

- At low angles, refraction occurs.
- Beyond the critical angle, the light reflects entirely within the block.

**Conclusion**

Total internal reflection occurs when the angle of incidence exceeds the critical angle. This confirms the principle of TIR.

---

## EXP.6

### Convex and Concave Lenses

**Objective**

To study the behavior of light rays through convex and concave lenses.

**Apparatus**

- Convex and concave lenses
- Ray box
- Screen
- Meter scale

**Procedure**

1. Shine parallel rays onto each lens and observe their paths.
2. Note focal point in convex (converging) and apparent divergence in concave (diverging).
3. Measure focal length by adjusting the distance until a clear image forms.

**Observations**

- Convex lens focuses rays at focal point.

- Concave lens diverges rays as if from a virtual focal point.

**Conclusion**

The experiment validates the converging nature of convex lenses and diverging nature of concave lenses. Their focal lengths can be measured through image formation.

---